

PW Analytical Balances



The PW series of Analytical Balances are designed to meet the needs of Laboratories, Schools, Industry and commercial users offering a wide range of applications.

Key to symbols used : Move mouse over a symbol for an explanation

Features

- Metal casing
- Force restoration mechanism
- RS-232 bi-directional interface
- Accessories interface
- Security locking point
- Internal calibration
- GLP print outs
- Date and time
- Colour coded keys
- Backlit LCD with dual display
- Up to 16 weighing units (including one custom unit)
- Capacity tracker
- Multilingual display

Applications

- Weighing
- Net / Total
- Check weighing
- Percentage weighing
- Weight accumulation
- Averaging
- Parts counting
- Animal / dynamic weighing
- Density determination
- Below balance weighing

Accessories

- No.8023 Adam Printer
- No.5001 Adam Statistics Printer
- No.3584 Thermal Paper Roll for Adam Printer
- No.8036 Anti-Vibration Table
- No.7452 Density Determination Kit
- No.3620 Below Balance Hanger
- No.8030 Security Lock and Cable
- No.8003 Dust Cover
- No.9028 RS-232 Cable
- No.9061 AdamDU - Data Collection Program

International Co. for Lab. Instruments

45 Shalhoub st. Ain Shams, Cairo, Egypt Post Code 11311



+202 24916139 or +202 2491613



+202 24953441