

## CPWplus Balance Bench Scales



The CPWplus is a multi-purpose industrial scale.

Key to symbols used : Move mouse over a symbol for an explanation

### Features

- Stainless steel weighing pan with remote indicator
- Large backlit LCD display
- Simple 4 button operation
- Portable with battery power
- Kg/ lb/ oz/ lb & oz switchable
- Tare function
- Mains adapter and wall mounting kit supplied as standard
- Hold function
- RS-232 bi-directional interface

### Applications

- Weighing
- Animal / dynamic weighing

### Accessories

- No.8023 Adam Printer
- No.7954 Hard Carry Case with Lock
- No.3584 Thermal Paper Roll for Adam Printer
- No.9014 RS-232 Cable
- No.9061 AdamDU - Data Collection Program

## CBK Bench Scales



Reliable performance, durable construction and versatility make the CBK unmatched for value in industrial weighing. Serving a variety of applications the CBK includes basic weighing, check weighing with LED limit indicators and parts counting. Full featured and easy-to-use, the CBK is the right choice for industry when flexibility matters.

Key to symbols used : Move mouse over symbols for an explanation

### Features

- EC type approval
- Large backlit LCD display
- Overload protection
- Easy to use keypad
- Rechargeable battery with up to 90 hours operation
- Tare function
- Large stainless steel pan
- Memory accumulation
- Pre-set weighing with alarm
- Check weighing LED bar graph
- Capacity tracker

### Applications

- Weighing
- Check weighing
- Percentage weighing
- Weight accumulation
- Parts counting

### Accessories

- No.8023 Adam Printer
- No.3584 Thermal Paper Roll for Adam Printer
- No.5006 In-use Cover
- No.5006-5 In-Use Cover pack of 5
- No.5006-10 In-Use Cover pack of 10
- No.6081 Weigh Below Hanger (factory fitted)
- No.9028 RS-232 Cable
- No.9061 AdamDU - Data Collection Program
- No.9025 Hard Carry Case with Lock

## QBW series Bench Scales



The QBW series is ideal for a wide range of weighing applications.

Key to symbols used : Move mouse over symbols for an explanation

### Features

- Large LCD display
- Large stainless steel pan
- Sealed keypad
- Rechargeable battery with up to 150 hours of operation / mains adapter supplied as standard
- Tare function
- g, kg, lb, oz and lb & oz
- Overload protection

### Accessories

- No.5006 In-use Cover
- No.5006-5 In-Use Cover pack of 5
- No.5006-10 In-Use Cover pack of 10
- No.9025 Hard Carry Case with Lock
- No.6081 Weigh Below Hanger (factory fitted)

---

**International Co. for Lab. Instruments**

45 Shalhoub st. Ain Shams, Cairo, Egypt+ Post Code 11311



+202 24916139 or +202 24916139



+202 24953441