



The PGW Series of Precision Balances are designed to meet the needs of Laboratories, Schools, Industry and commercial users offering a wide range of applications.

Key to symbols used : Move mouse over a symbol for an explanation

Features

- Metal casing
- Force restoration mechanism
- RS-232 bi-directional interface
- Accessories interface
- Security locking station
- External calibration
- Internal calibration models
- GLP print outs
- Date and time
- Colour coded keys
- Backlit LCD with dual display
- Up to 18 weighing units (including one custom unit)
- Draft shield supplied as standard with models with 0.001g readability
- Capacity tracker
- Multilingual display

Applications

- Weighing
- Net / Total
- Check weighing
- Percentage weighing
- Weight accumulation
- Averaging
- Parts counting
- Animal / dynamic weighing
- Density determination
- Below balance weighing

Accessories

- No.8023 Adam Printer
- No.5001 Adam Statistics Printer
- No.3584 Thermal Paper Roll for Adam Printer

- No.7453 Density Determination Kit (mg models only)
- No.8036 Anti-Vibration Table
- No.3620 Below Balance Hanger
- No.8030 Security Lock and Cable
- No.8002 Dust Cover
- No.8031 In-use Cover
- No.9028 RS-232 Cable
- No.9061 AdamDU - Data Collection Program

ADK series Precision Balances



The ADK series are high capacity and high accuracy precision balances which are ideal for many applications such as parts counting, percentage weighing, animal weighing, etc.

Key to symbols used : Move mouse over symbols for an explanation

Features

- 3 weighing units - g, lb & oz
- RS-232 bi-directional interface
- Full range tare
- Weigh below facility with optional hanger
- Robust construction with heavy duty casing
- Large stainless steel weighing pan
- Vacuum Fluorescent Display
- Simple key operation
- Anti-shock protection

Applications

- Weighing
- Percentage weighing
- Below balance weighing
- Parts counting
- Animal / dynamic weighing

Accessories

- No.8023 Adam Printer
- No.3584 Thermal Paper Roll for Adam Printer
- No.8036 Anti-Vibration Table

- No.3620 Below Balance Hanger
- No.9028 RS-232 Cable
- No.9061 AdamDU - Data Collection Program

International Co. for Lab. Instruments

45 Shalhoub st. Ain Shams, Cairo, Egypt Post Code 11311



+202 24916139 or +202 24916139



+202 24953441