

ACBplus compact balance



The ACBplus is a lightweight, rugged precision compact balance suitable for a wide range of weighing applications.

Key to symbols used : Move mouse over symbols for an explanation

Features

- Capacity tracker
- Rechargeable battery / mains adapter supplied as standard
- RS-232 bi-directional interface
- Large backlit LCD display
- 15 weighing units
- Auto power off
- Stainless steel pan
- Draft shield supplied as standard
- Integral security bracket
- Adjustable feet

Applications

- Weighing
- Percentage weighing
- Weight accumulation
- Parts counting
- Below balance weighing

Accessories

- No.8023 Adam Printer
- No.3584 Thermal Paper Roll for Adam Printer
- No.5111 In-use Cover
- No.7958 Hard Carry Case with Lock
- No.9028 RS-232 Cable
- No.9061 AdamDU - Data Collection Program

AQT compact balance



The AQT is an economically priced balance tailored for many weighing applications.

Key to symbols used : Move mouse over symbols for an explanation

Features

- Stainless steel pan
- Clear backlit LCD display
- Auto power off
- Mains or battery operated
- Mains adapter supplied as standard
- Several weighing units
- Draft shield supplied as standard
- Integral security bracket
- Adjustable feet

Applications

- Weighing

Accessories

- No.5111 In-use Cover
- No.7958 Hard Carry Case with Lock

International Co. for Lab. Instruments

45 Shalhoub st. Ain Shams, Cairo, Egypt+ Post Code 11311



+202 24916139 or +202 24916139



+202 24953441