

Laminar Flow Cabinets

AURA HZ

CROSS-FLOW CABINETS: Quality with a flair of tradition

Quality with a flair of tradition

- Horizontal laminar air flow cabinet.
- Last generation Microprocessor control
- Automatic airflow monitoring
- Reduced air flow switch
- H14 filter with micromesh downstream equalising plenum for the highest airflow speed uniformity
- Micromesh foil for maximum eye comfort.
- Side glass walls with Lateral Flow Effect (LFCTM) assembly.
- Silent and quite operation <59dB(A)

Available in two sizes: 1.2 and 1.8 mt



AURA VERTICAL SD4

DOWN-FLOW CABINETS: Quality with a flair of tradition

- Vertical laminar air flow cabinet.
- Sliding Sash front window
- Electronic Control
- H14 filter with micromesh downstream equalising plenum for the highest airflow speed uniformity
- Side glass walls.

- Soft touch keys



AURA MINI

Downflow CABINETS: When space is at a premium

Compact, vertical laminar air flow cabinet.

- Reduced size for the easiest installation in small or crowded labs.
- Easily convertible inward or outward barrier operating modes.
- Hinged front panel and side panels in tempered glass
- H14 downflow filter
- Filtrete® exhaust filter (or prefilters depending on configurations)
- Active PCR version
- Silent and quite operation <60dB(A)



Code	Description	Format	Brand	
LH20000	AURA HZ 48 Horizontal Laminar Flow Cabinet, 120 cm wide.	1 pz	BIOAIR	
LV30000	AURA Mini - Vertical Laminar Flow Cabinet, cm 75 wide.	1 pz	BIOAIR	
LV40000	AURA VERTICAL SD4 Vertical Laminar Flow Cabinet, cm 124 wide	1 pz	BIOAIR	

International Co. for Lab. Instruments

45 Shalhoub st. Ain Shams, Cairo, Egypt Post Code 11311



+202 24916139 or +202 24916139



+202 24953441