

Machine Specifications for Model CP-501 Single Tablet Press

Accessories Included

- Diameter 8mm Convex Round shaped punches and die installed
- Tool box
- Sample of tablets
- Sample of placebo powder is included

Specifications

- cGMP Standard Construction
- The machine is completely set up and ready for operation upon arrival
- Variable speed control
- Overload protection function
- Safety Interlock Systems
- It makes round shaped, other irrigable shaped and even ring-like tablet
- Easy changing parts and cleaning
- Ergonomic and Rugged Design
- Maximum Pressure: 30 kn
- Maximum Diameter of Tablet: 20mm
- Maximum Depth of Fill: 16mm
- Maximum Thickness of Tablet: 8mm
- Best single tablet press in the industry

Capacity

- Up to 10,800 tablets per hour

Package Dimensions & Weight

- Machine Dimensions: 34"x 26" x 35"
- Net Weight: 330 pounds
- Gross Weight: 435 pounds
- Package: 01 Crated Wooden Box
- Packing Dimensions: 47" x 28" x 40"

Power Source

- A.C. 220V, 50Hz, Single Phase

Warranty

Valid 12 months from date of shipment



International Co. for Lab. Instruments

45 Shalhoub st. Ain Shams, Cairo, Egypt Post Code 11311



+202 24916139 or +202 24916139



+202 24953441